## Hot/Cold Pack

### INSTRUCTIONS FOR USE:

- Provides effective relief from muscular pain and helps reduce swelling.
- Great for treating a wide range of soft tissue injuries.
- Pack will remain flexible at low temperatures.
- Pack lasts longer and is more effective than a standard ice or hot pack.
- Reusable pack is perfect for continuous treatment at home.

#### For cold use:

Place pack in the freezer for up to 2 hours before use. Wrap in a cloth or protective sleeve before applying directly to the skin. Use pack in short intervals no longer than 20 mins.

#### For hot use:

Place pack in the microwave for 20 seconds at 600W. Increase heating time in 10 second intervals to reach desired temperature if needed. Times should be reduced for microwaves with higher power ratings. Wrap in a cloth or protective sleeve before applying directly to the skin. Use pack in short intervals no longer than 20 mins.

#### Warning:

Never place pack directly onto the skin as extreme temperatures can cause burns. Always wrap pack in a cloth or protective sleeve before use. Always use under adult supervision. Keep out of the reach of children. If you are affected by circulatory problems seek medical advice before use. Cease use if pack causes any discomfort. If the pack is punctured or opened and the solution is accidentally swallowed, drink plenty of water. If solution contacts eyes, wash out thoroughly with water. Consult your doctor if you have any concerns about your health. Check pack for leakage before use. Clean pack before use. This pack contains Sodium Carboxymethyl Cellulose, Glycerol and Water.

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Henry Schein U.K. Holdings Ltd.